


Install hCare Portal on your Android device.

This section details how to install the Handheld Native application on an Android device running Android OS 2.3.4. Some screens may be different based on the OS version.

Install the Handheld Native Platform


1. Press the applications button on the phone. Navigate to Settings. Select Applications. Check Unknown Sources.




Menu>Settings>Applications

In the address field, type the name of your mobilizer server appended with <https://dedv.mobile.medcity.net/ota>


and press [Enter].




2. Once the download is complete, press the Home button on the device. From the home screen pull down the Android Notification screen.




3. Click on the PatientKeeper.apk file.



4. Click the Install button.



5. After the application is installed click Done.



6. To add the icon to the home screen press the application button from the home screen, find the PatientKeeper icons, press and hold the PatientKeeper icon, and drag it to a place on the screen. This will create a shortcut to the application.



Configure the Application

Follow the steps below to configure the application when you use it for the first time.

1. Click the PatientKeeper icon you added in the previous procedure to launch the application.



2. In the Host field enter the division-specific mobilizer

<https://dedv.mobile.medcity.net> .

Host:

Use Storage Card

MOBILIZER_7_4_1-20120125_1524

Test	Logs
Save	Cancel

3. Press the Save button.

Host:

Use Storage Card

MOBILIZER_7_4_1-20120125_1524


Test	Logs
Save	Cancel

4. Enter the Username and Password. Press Login.


Login to Server

Username	<input type="text"/>
Password	<input type="password"/>

Login



5. Device will sync. Patient List is displayed.



No Profile Selected	
UPDATED: 3/1 8:44 am	
All Patients (25)	
ADM, INPT5	26Y M
J.CAMPUS4.J.CAMPUS4.8.COCQA1A	
TEST	
ALBRIGHT, ADMCRIT	F
MR.1N.MR.10.1.COCMA1A	
TEST PHAR ADM CRITIERIA INFO	
ALBRIGHT, COCMA1A	43Y F
MR.1E.ROOM2DEC20.BD1.COCMA1A	
TEST	
ALBRIGHT, COMPLETE	42Y F
MT.8E.MT8.1.COCMA2	
MR# MT00000059-COCMA2-TEST M...	
ALBRIGHT, DEFECT160	42Y F
MR.1E.ROOM DEC11.BD1.COCMA1A	
TEST ADMIT OF CANCELLED ACCT D...	
ALBRIGHT, OBS	42Y F
MR.1N.MR.10.11.COCMA1A	
TEST OBS PATIENT	